3C-BOND

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3C-BOND: micro-fine (0,01 to 0,02 mm) bonding agent to fire porcelain to constructions of conventional chrome-cobalt alloys

Field of application: crown and bridgework and removable Cr-Co constructions

3C-BOND prevents the yellow-green discolouration of porcelain, especially in the thin margins with Ni-Cr alloys.

INSTRUCTIONS FOR USE

Metal preparation

Sandblast with aluminium oxide. Clean in an ultrasonic bath or by boiling for 10 minutes in distilled water, or steam clean.

The metal frame must be completely dry and should only be touched with clean tweezers.

Mix 3C-BOND with distilled water and apply in a thin layer. Cover also marginal areas.

Never re-use already mixed 3C-BOND.

Fire 3C-BOND under vacuum in the same way as opaque.

Preheating temperature: 650 °C (1200 °F)
Temperature rise: 55 °C (131 °F) per min
Final temperature: 980 °C (1800 °F)

After firing, immediately remove the metal frame from the furnace and bench cool.

3C-BOND should now have a yellowish, slightly shiny surface.

Continue with the opaque firing cycle as recommended by the porcelain manufacturer.

TROUBLESHOOTING

Before use

Before firing 3C-BOND to the metal frame, we recommend a firing test on a Cr-Co ingot or casting cone to examine the reaction of 3C-BOND with the Cr-Co alloy to be used.

Green-black stains

The 3C-BOND mixture was too thin.
- The metal frame was still greasy or wet.
- 3C-BOND was mixed more than once.

In these cases, thoroughly clean the metal frame and start again.

3C-BOND shows bubbles.

Degas the metal frame under vacuum for 10 minutes prior to repeated application of 3C-BOND.

3C-BOND shows cracks.

The 3C-BOND mixture was too thick. Once again blast and thoroughly clean. Then apply and fire 3C-BOND.

Cracks in the porcelain

If 3C-BOND is rust-coloured instead of yellow after firing the bonding agent, the Cr-Co alloy that was used is not suitable or 3C-BOND was re-used.

SAFETY INSTRUCTIONS

Skin contact

Immediately wash with water and soap. 3C-BOND normally does not irritate the skin.

Eye contact

Carefully flush the eyes with running water.

Ingestion

Seek medical advice if irritation persists.

STORAGE AND SHELF LIFE

3C-BOND must be stored in a well-sealed package in a dry place.

In case of proper storage, the product has an almost unlimited shelf life.

The package may be disposed of as household waste.

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Copyright© ALPHADENT NV 8/2016