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ADVANTAGES / FEATURES

CEKA SITE is a bonding and filling material that is mainly used in the attachment technique. The anaerobic properties of CEKA SITE allow to bond components whilst expelled material sets slower. Outstanding bond values to precious and non-precious metals. No light curing necessary. Simple processing by sandblasting first. Excellent setting behaviour. Easy to apply. Highly flowable material. Easy clean-up.

INSTRUCTIONS FOR USE

Always use CEKA SITE at room temperature. Sandblast the metal surfaces to be bonded with coarse aluminium oxide. Use the automix tip to expel CEKA SITE. You now have maximum 1.5 minutes working time. CEKA SITE is sufficiently set after 10 minutes. The expelled excess sets slowly as exposed to air and thus cannot serve as a control. Be sure to make a new mixture of CEKA SITE for each use. Use optional CEKA Multi Primer for even higher bond values.

SIDE-EFFECTS

Direct contact with uncured CEKA SITE may cause skin irritation and allergic reactions.

STORAGE AND SHELF LIFE

Limited shelf life: see date on label.
Store in a cool place between 4 °C and 10 °C (39 °F and 50 °F).

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