

PRECI-BALL

A MULTIPLE USE ATTACHMENT

2.25 mm
ball diameter



A 2.25 mm ball diameter makes PRECI-BALL compatible with most implant systems. Combine with a ball abutment of the implant company.



The ball may be incorporated into a bar.



The PRECI-BALL female is also available in a bulk package of 20 for great savings.



The PRECI-BALL system is fully interchangeable with PRECI-CLIX (see PRECI-CLIX leaflet).

TABLE OF ALLOYS

ORAX - Yellow - Au 70 - Pt 5 - Ag 12 - Cu 13

Melting range: 910-955 °C (1670-1823 °F)

PALLAX - White - Au 2 - Ag 37 - Pt 9.5 - Pd 37 - Cu 12.5 - Co 2

Melting range: 1055-1130 °C (1931-2066 °F)

TITANAX - White - Ti 90 - Al 6 - V 4 (Do not heat!)

IRAX - White - Au 60 - Pt 24 - Pd 15 - Ir 1

Melting range: 1400-1460 °C (2552-2660 °F)

WWW.CEKA-PRECLINE.COM

SWISS MADE



CEKA REVAX
CEKA SOL
PRECI-HORIX
PRECI-VERTEX
PRECI-VERTEX AT
PRECI-BAR
PRECI-BALL
PRECI-CLIP
PRECI-SAGIX
PRECI-CLIX AXIAL
PRECI-CLIX EXTRACORONAL
PRECI-CLIX RADICULAR
PRECI-CLIX RADICULAR RC
PRECI-TUBIX/MORTIX
PRECI-POST
PRECI-PROFILE
PERMA-RET
3C-BOND
EXPANDO
CEKA-MAX

PRECI-BALL

PRECI-BALL

CEKA
ATTACHMENTS
PRECI-LINE

Prospectus no 267.E-ed.C
Copyright© ALPHADENT NV 11/2009
BTW BE 0414.789.321 RPR Kortrijk

ARSEUS DENTAL
SERVING THE DENTAL PROFESSIONALS

ALPHADENT NV

Textielstraat 24

8790 Waregem, Belgium

tel +32 56 629 100

fax +32 56 629 101


e-mail: info@ceka-preciline.com



PRECI-BALL

AN ATTACHMENT WITH OPTIONS


1 Solder the male part (gold or palladium alloy) into the post-coping. This combines true precision and accuracy of technique at a low cost.



2 Cast-on with a precious or palladium-based alloy and have the male part replaceable. Recommended technique to combine precision and serviceability.



3 Cast in a hard dental alloy. The most economical solution. Use the PRECI-BALL cup burr to finish the ball.



Relieve the post-coping with the supplied tinfoil space maintainer.
1201 E

PRECI-BALL STANDARD AU # 1201

Female in ORAX, for retention in acrylic resin

1201 A 1 female
1201 AB 20 females

Male in ORAX, for soldering to the post-coping

1201 C 1 male

H: 3.5 mm

PRECI-BALL STANDARD TI # 1203

Female in TITANAX, for retention in acrylic resin

1205 A 1 female
1205 AB 20 females

Male in PALLAX, for soldering to the post-coping

1203 C 1 male

H: 3.5 mm

PRECI-BALL PC # 1205

Female in TITANAX, for retention in acrylic resin

1205 A 1 female
1205 AB 20 females

Male in burn-out plastic, for casting with the post-coping

1205 C 1 male
1205 CB 50 males

H: 3.5 mm

PRECI-BALL CV # 1206

Female in ORAX, for retention in acrylic resin

1201 A 1 female
1201 AB 20 females

Male in TITANAX, to thread into the base ring

1206 C 1 male

Base ring in IRAX, for cast-on with precious or palladium-based alloys

RA 0063 1 base ring

H: 4.3 mm

ACCESSORIES

	Activating tool # 1220	1 piece
	Desactivating tool # 1221	1 piece
	Screw-driver for the threaded male # IMP-XS-042	1 piece
	Paralleling mandrel for males # 1201 P	1 piece
	Paralleling mandrel for females # 1211 P	1 piece
	Paralleling mandrel for base rings # RE P 4	1 piece
	Cup burr # 1230	2 pieces Burr to finish the 1205 C male
	Model analogue # 1201 D	2 pieces
	Male V extra long # 1206 CL	1 piece Special titanium male with 0.3 mm longer threading
	CEKA BOND # CB1	5 ml Anaerobic adhesive to prevent gradual unthreading of prosthetic components