

PROCESSING PERMA-RET

1. Determine the exact position for the **1602** PERMA-RET wire.
2. Place the **1606** burr into a handpiece. Drill an approximately 3 mm deep hole in the tooth.
3. Place the retention wire in the **1604** wire holder.
4. Screw the wire into the opening.
5. Cut the wire to the required length (approximately 4 mm) with pliers or a separating disc.
6. Bend or squeeze the end of the wire. It will then offer retention for acrylic resin.

Side effects, warnings and precautions

- The attachments are intended for single use.
- The products are non-sterile.
- There is a risk of poor fit when patient conditions change.
- Bacterial adhesion can be avoided by applying hygiene measures.
- Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
- The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
- For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
- Do not heat items containing titanium.
- Do not use items containing nickel in case of nickel allergy.
- The accessories RE H 79 and H 35 must be used outside the mouth.