

## **BONDING A COMBINATION WORK**

The working procedure and the precision for bonding the abutment retainers are extremely important. They guarantee the parallelism and a longer life of the combination work.

1. Check if the construction parts fit accurately in the mouth.
2. Lubricate the space maintainers with vaseline and place them on the spring pins.
3. Bond the abutment retainers.
4. Immediately press the removable prosthesis over the abutment retainers and bring the patient's jaws manually in centric relation.
5. Allow the composite to set.
6. Remove the secondary construction and the composite residues from the primary and secondary parts.
7. Check the snap and the retention of the attachments (see INFO 059 and INFO 069).