

CONVERSION OF PRECI-VERTIX AT INTO PRECI-VERTIX

1. Fill up the impression tray with impression material and take an impression.
2. Place the male analogue **1819** in the impression tray with the retention opening facing up.
3. Pour the model in stone.
4. Pour the plastic housing **1808** in an alloy of your choice.
5. Remove the entire female from the prosthesis. It is cemented and can therefore easily be removed with a burr. Provide sufficient space to incorporate the new housing **1808**.
6. Reduce the cast housing, if necessary. For improved adhesion, additional retention can be created on the extension of the housing.
7. You now have the option to bond the cast housing with **CEKA SITE** or to polymerize it in acrylic resin.
8. Position the female **1802**, **1805** or **1806** in the housing.

Side effects, warnings and precautions

- The attachments are intended for single use.
- The products are non-sterile.
- There is a risk of poor fit when patient conditions change.
- Bacterial adhesion can be avoided by applying hygiene measures.
- Inappropriate use or bad manufacturing can lead to premature wear of the attachments.
- The functionality of the attachments will be adversely affected by traumas such as grinding and bruxism.
- For the purpose of traceability we advise you to record the lot number of the applied products in the patient file.
- Do not heat items containing titanium.
- Do not use items containing nickel in case of nickel allergy.
- The accessories RE H 79 and H 35 must be used outside the mouth.