

# SPECIALTIES

## CEKA SOL ENGLISH



CEKA  
ATTACHMENTS  
PRECI-LINE

# CEKA SOL

FLUX INCORPORATED



# LOW-FUSING TUBE SOLDER WITH FLUX INCORPORATED

A complete range of low-fusing solders for multiple applications. The incorporated flux provides for easy soldering of dissimilar metals.



*Simple & easy!*

PRECIOUS  
NON-PRECIOUS  
CR-CO

**CEKA SOL # CS00** X X X

Yellow 780-820 °C (1436-1508 °F)  
Au 80 - Cu 13 - In 5 - Zn 2

**CEKA SOL FILIGRAN # CS00 F** X X X

Yellow 780-820 °C (1436-1508 °F) (contains less flux)  
Au 80 - Cu 13 - In 5 - Zn 2

**CEKA SOL W # CS22** X X X

White 865-930 °C (1589-1706 °F)  
Au 76.5 - Cu 0.4 - Zn 6.5 - Ni 16.6

**CEKA SOL L # CS21** X X

Yellow 615-630 °C (1139-1166 °F)  
Au 42 - Ag 38.5 - Zn 11 - In 8.5

**CEKA SOL G # CS20** X

Yellow 755-820 °C (1391-1508 °F)  
Au 68 - Ag 13.3 - Cu 12.8 - Sn 1 - Zn 4.9

**CEKA SOL CR # CS23** X

White 1000-1030 °C (1832-1886 °F)  
Au 10 - Ag 37 - Pd 37 - Cu 12 - Zn 4

## ADDITION ALLOY FOR LASER WELDING WITHOUT FLUX

**CEKA SOL T1 # CS31**

White, pure titanium  
Round 8 x 25 cm Diameter: 0.4 mm

**CEKA SOL OR # CS32**

Yellow 910-995 °C (1670-1823 °F)  
Au 70 - Pt 5 - Ag 12 - Cu 13  
Round 2 x 25 cm Diameter: 0.4 mm

**CEKA SOL PA # CS33**

White 1055-1130 °C (1931-2066 °F)  
Au 2 - Pt 9.5 - Ag 37 - Pd 37 - Cu 12.5 - Co 2  
Round 2 x 25 cm Diameter: 0.4 mm



CEKA AXIAL  
PRECI-CLIX AXIAL  
PRECI-CLIX RADICULAR  
PRECI-BALL

**01.**  
AXIAL  
ATTACHMENTS

CEKA EXTRACORONAL  
CEKA ID  
PRECI-VERTIX  
PRECI-VERTIX AT  
PRECI-CLIX EC  
PRECI-SAGIX  
PRECI-52  
PRECI-TUBIX/MORTIX

**02.**  
EXTRACORONAL  
ATTACHMENTS

PRECI-BAR  
PRECI-HORIX  
PRECI-CLIP

**03.**  
PROFILE  
ATTACHMENTS

PRECI-PROFILE  
PRECI-POST  
**CEKA SOL**  
CEKALLOY  
CEKA SITE  
CEKA BOND  
3C-BOND  
PERMA-RET  
PRECI-SEP  
PLASTICWAX  
EXPANDO

**SPECIALTIES**

[WWW.CKPL.EU](http://WWW.CKPL.EU)

CEKA  
ATTACHMENTS  
PRECI-LINE



Prospectus n° 268.E/D/N-ed.A

© Copyright ALPHADENT NV 10/2011 - BTW BE 0414.789.321 RPR Kortrijk

ALPHADENT NV • Textielstraat 24 8790 Waregem, Belgium

T +32 56 629 100 - F +32 56 629 101 - E info@ckpl.eu

