


**PRECI-HORIX STANDARD**  
**# 1701**  
 2 plastic bar patterns  
 6 yellow riders  
 6 processing riders  
 1 insertion tool

**KIT**



**PRECI-HORIX PLASTICWAX**  
**# 1731**  
 2 PLASTICWAX bar patterns  
 6 yellow riders  
 6 processing riders  
 1 insertion tool

**KIT**



**PRECI-HORIX COMBI**  
**# 1741**  
 1 plastic bar pattern  
 2 yellow riders  
 2 processing riders  
 2 stainless steel housings  
 1 insertion tool

**KIT**



**PRECI-HORIX RIDER**  
 Normal retention:

**# 1802 yellow 6 pcs**  
**# 1802B yellow 50 pcs**

Reduced retention:

**# 1805 white 6 pcs**  
**# 1805B white 50 pcs**

Increased retention:

**# 1806 red 6 pcs**  
**# 1806B red 50 pcs**

Combination package:

**# 1807 2 y + 2 r + 2 w**  
**# 1807B 25 y + 15 r + 10 w**



**PRECI-HORIX METAL HOUSING**

Ultrathin stainless steel housing with titanium/nitrate layer (gold-coloured) for a perfect fit of the riders in the prosthesis.


**# 1706 6 pcs**  
**# 1706B 50 pcs**



**PRECI-HORIX PLASTIC PATTERN**

Burn-out plastic pattern.  
 Cast in any hard alloy.  
 Length: 50 mm - Height: 4 mm - ø 1.8 mm

**# 1704 6 pcs**  
**# 1704B 50 pcs**



**PRECI-HORIX PLASTICWAX PATTERN**

Ideal for use with fast set investments.  
 Cast in any hard alloy.  
 Length: 50 mm - Height: 8.25 mm - ø 1.8 mm

**# 1734 6 pcs**  
**# 1734B 50 pcs**

CEKA AXIAL  
 PRECI-CLIX AXIAL  
 PRECI-CLIX RADICULAR  
 PRECI-BALL

CEKA EXTRACORONAL  
 PRECI-VERTIX  
 PRECI-CLIX EC  
 PRECI-SAGIX  
 PRECI-52

PRECI-BAR  
**PRECI-HORIX**  
 PRECI-CLIP

PRECI-PROFILE  
 PRECI-POST  
 CEKA SOL  
 CEKA SITE  
 CEKA BOND  
 3C-BOND  
 PERMA-RET  
 PRECI-SEP  
 PLASTICWAX  
 EXPANDO  
 CEKA Multi

**01.**  
 AXIAL  
 ATTACHMENTS

**02.**  
 EXTRACORONAL  
 ATTACHMENTS

**03.**  
 PROFILE  
 ATTACHMENTS

SPECIALTIES

**03.**  
 PROFILE  
 ATTACHMENTS

**PRECI-HORIX®**  
 ENGLISH

[WWW.CKPL.EU](http://WWW.CKPL.EU)



**CEKA PRECI-LINE APP**



Prospectus n° 273.F.ed.C  
 © Copyright ALPHADENT NV 9/2014 · BTW BE 0414.789.321 RPR Kortrijk  
**ALPHADENT NV** · Mannebeekstraat 33 8790 Waregem, Belgium  
 T+32 56 629 531 · E info@ckpl.eu


**CEKA**  
 ATTACHMENTS  
**PRECI-LINE**

**CEKA**  
 ATTACHMENTS  
**PRECI-LINE**

# PRECI-HORIX

SIMPLE AND ECONOMICAL



The PRECI-HORIX bar and rider system, developed by Helmut Hader, is the most popular of all bar systems due to its economy and simplicity. The excellent retention of the plastic clips, the patient comfort, the ease of processing and durability make PRECI-HORIX a first choice by experts for over 40 years.



- ✓ Easily replaceable female riders in three retention levels to choose from
- ✓ The bar pattern may be cast in alloy of choice
- ✓ Extremely suitable for use on implants
- ✓ Direct support and rotational movement
- ✓ The soft snap of the female riders offers the patient security



## PRECI-HORIX HOUSING SIMPLIFIES CLIP PLACEMENT



Recommended use is processing the metal housings into the denture simplifying clip placement and then using one of the three retention grips of plastic PRECI-HORIX clips in the metal housings.

White: reduced retention  
Yellow: normal retention  
Red: increased retention



The prevailing 1.8 mm bar pattern is also compatible with PRECI-CLIP riders.

## EASY PROCESSING

